

Innovation and Creativity

Technical support and Equipment in Process Metallurgy

Delphine Leroy – Senior engineer IMT SiManTi 3.10.2013

Technical support and equipment in Process Metallurgy

- Technical support at IMT
- Technical services at NT faculty
- Equipment in process metallurgy
 - Lab scale, small scale, pilot
 - NTNU, Sintef
- Procedure when running experiment
- Wiki database for equipment



Technical support at IMT

- 15 technical employee at IMT
- 7 in Alfred Getz vei (Agv)
- 8 in Kjemi blokk II (KII)
- <u>Technical staff IMT</u>
 http://www.ntnu.edu/mse/employees/technical-staff



Engineer in process metallurgy

- High temperature equipment : graphite tube furnaces, induction furnaces...
- Crushing/sieving lab
- Instrumentation: thermocouple (K,S,B and C), MFC, datalogger...
- Diverse lab consumables : crucibles, tube, hoses...
- Room responsible : HSE



Electron microscope laboratories

Senior Engineer Morten P. Raanes

- Electron probe micro-analyser EPMA
- EM Lab
- http://www.material.ntnu.no/lab/material/booking.php



Senior Engineer Yingda Yu

- Scanning electron microscope
 SEM
- Transmission electron microscopy TEM
- Liquid nitrogen





Staff engineer Torild Krogstad

- XRD
- Etching lab
- GD-MS, Glow discharge mass spectrocropy
- ICP-MS, inductive coupled plasma
- Chemical, <u>Eco Online</u>
 https://innsida.ntnu.no/wiki/-/wiki/Norsk/Stoffkartotek
- Gas
- Metallography lab





Senior Engineer Trygve Schanche

- Metallography lab Agv and KII
- Sample preparation
- Heat treatment
- Cutting
- Microscope
- Hardness measurement



Staff engineer Pål Chr. Skaret

- Mechanic testing
- Hydrolic press
- Rolling and torsion
- Officices
- PC
- Telephon
- Access card
 - Risk assement
 - Room access requisition
 - K:\imt\HMS





Staff engineer Harald Holm

- Electric and electronic support
- Safety representative Agv
- https://innsida.ntnu.no/wiki/-/wiki/English/Safety+represent ative
- PC
- Keys



Staff engineer Gunn Torill Wikdahl

- HSE coordinator
- New PhD ? Get HSE course
- HMS at NT faculty
- https://innsida.ntnu.no/wiki/-/wiki/Norsk/HMS+ved+NT
- K:\imt\HMS





Staff engineer Solveig Jonassen

- Simple purchases
 - Supplier agreement
 - E-requisition
 - https://innsida.ntnu.no/bestille
- Chemical storage, EcoOnline KII
- Room access





- Senior Engineer Eli Beate Larsen
 - Technical staff leader IMT, Thermic analysis, Cold isostatic press



- Senior Engineer Elin H. Albertsen
 - Carbon lab, hazard waste



electrochemical laboratories





- Staff Engineer Pei Na Kui
 - Furnaces KII, glove box, spray pyrolyse lab



- Gas central system KII, Light microscope



- Senior Engineer Julian Tolchard
 - XRD and SEM in KII



Technical services at the NT faculty

- Glass Blowing Workshop
 - Realfagbygget DU3
- Mechanical/Engineering Workshop
 - Realfagbygget DU3
- IT and Instrumentation Service
 - Realfagbygget C2
- Order work from innsida / mine bokmerker / NT Fak.
- Order work at NT workshop http://www.nt.ntnu.no/adm/Verksted/
- <u>Description of the workshops</u>
 http://www.ntnu.no/nt/om-oss/verksteder



Glass Blowing Workshop

Manager: Astrid Salvesen



- Glass equipment manufacturing and repair. Special design.
- Various qualities of glass and materials:
 - Borosilicate glass (Duran glass)
 - Quarts (SiO₂)
 - Ceramic materials
 - Noble metals (gold, platinum and various alloys)
- Machining: diamond sawing, drilling, grinding and polishing



Mechanical/Engineering Workshop Manager: Harald Snekvik



- Construction, development, modification, repair and maintenance of apparatus and equipment using machining techniques
 - Machining of various materials (high alloy qualities of steel, plastic materials, wooden materials, etc).
 - Graphite (Aksel Alstad), coke, boron-nitride and various composite materials
 - Modern CNC turning lathes and milling machines. EDM equipment (spark erosion). Welding (traditional and plasma), soldering, etc
- Competence on high vacuum technology
- Technical calculations and drawings (AutoCad) according ISO and DIN standards



Innovation and Creativity

Equipment for metallurgy and mineral research

Type of equipment

- Characterisation
- Standardised tests
- Laboratory tests
- Small scale tests
- Pilot scale tests
- Industrial measurements

Types of research

- Mineralogy
- Ore dressing
- Powder technology
- Ferroalloy production
- PV Si production
- PV Si characterisation
- Alloy characterisation



Graphite tube furnaces

- Inert or reducing atmosphere
- Up to 2500°C
- Various features :
 - Weighting device
 - Gas lance
 - Quenching chamber
 - Translation of crucible
 - Off gas measurement
- Crucible size
 - Up to Ø45od*70
 - Resina Ø60/52*100







Vacuum induction furnace "Blue furnace"

- Temperature limit: 2000°C
- Manual heating control
- Crucible size
 - OD Ø80mm (normally)
 - Height 160mm, max 200mm
- Gas purging
- Sampling
- Crucible visible through glass






Cold Crucible induction furnace

- 75 kVA high frequency power generator (750 kHz)
- Multi wavelength pyrometer FAR FMP2
- Evacuating system with possibility of purging gas
- Freeze lining
- Work possibilities with very corrosive slags
- 300-500g






IF75 - Induction furnace

- Small scale smelting apparatus
- Investigation of reactant behaviour at different temperatures
- Volume of approx. 6 dm3
- Graphite crucible
- 75kW
- 3kHz





IF75 experiment example







Sample preparation lab


- Crushing
- Sieving
- Ultrasound bath
- Drying cabinet
- Pelletizing
- Mixing






DisVarDri furnace


- Sintef Shawn Wilson, Edith Thomassen
- Simulation of DRI shaft.
- Controlled gas inlet and temperature
- Online weight and gas composition measurement (IR)
- Tmax = 1200°C
- Sample 300g



Wettability furnace


- Sintef Tone Anzjøn
- Measurement of the contact angle of a liquid drop on a 10mm diameter substrate
 - The wetting angle
 - The melting point
 - The reactivity between different materials
- Heating rate up to 1600°C in 1 minute
- Atmosphere : vacuum, inert or reducing
- Maximum temperature 2400°C







Melting behaviour at increasing temperature




Pilot scale Si-production in 160kW single phase furnace


- Sintef Steinar Prytz
- 440 kVA power supply
- Possible to operate at different frequencies:
 - 25, 50, 75, or 100 Hz AC
- Automatic logging of power supply parameters and gas temperature.
- Opportunities for :
 - Off gas analysis (CO, CO₂, O₂, NO, NO₂ and SO₂)
 - Material balance
 - Energy balance





Pilot furnace Post experimental investigations





Samples for Microprobe, SEM and chemical analysis



How do I run experiment?

- Agree with your supervisor and write the risk assessment
- Get room access. Pål Skaret and room responsible
- Purchase <u>through the ordering system</u>:
 - Consumables. Delphine Leroy. Remember long delivering time!
 - Register chemical. Torild Krogstad. No storage in the office!
 - Gas. Torild Krogstad. 2 days delivering time.
- Working hours: 8:00 to 16:00, Monday to Friday
 - Contact supervisor and room responsible for late work
- Indicate what your doing on lab board
- Mark your samples
 - Tidy after yourself



Important documents in the lab

Safety card (vernekort)	Bills at entrances showing names of the Head of dept, HSE-coordinator, safety representative
Room card (romkort)	Bills outside entrances to laboratories showing vital information on safety issues to fire fighters and others
Apparatus card (apparaturkort)	Bills situated by the instruments telling who is responsible, what to do in emergency situations, etc.
Operating instructions (driftsinstruks)	How to perform a working task or using an instrument. Copies available in the lab.
Risk assessments (risikovurderinger)	Documentation of risks. Copies available in the lab.
Safety data sheets (sikkerhetsdatablad)	Safety information on dangerous substances. Copies available in the lab.



Good luck with your experiments!

